



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

CD234 Antibody [#G081]

Catalog Number: G081

Concentration: 1mg/ml

Swiss-Prot No.: Q16570

Other Names: Duffy antigen/chemokine receptor; Fy glycoprotein; GpFy; Glycoprotein D; Plasmodium vivax receptor; CD234 antigen; DARC; FY, GPD

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CD234.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CD234 antibody detects endogenous levels of total CD234 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:40000

References:

Chaudhuri A., Proc. Natl. Acad. Sci. U.S.A. 90:10793-10797(1993).

Tournamille C., Nat. Genet. 10:224-228(1995).

Iwamoto S., Blood 85:622-626(1995).

