



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

FSHR Antibody [#G101]

Catalog Number: G101

Concentration: 1mg/ml

Swiss-Prot No.: P23945

Other Names: Follicle-stimulating hormone receptor; FSH-R; Follitropin receptor; FSHR; LGR1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human FSHR.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FSHR antibody detects endogenous levels of total FSHR protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:1000

References:

Minegish T., Biochem. Biophys. Res. Commun. 175:1125-1130(1991).

Tilly J.L., Endocrinology 131:799-806(1992).

Kelton C.A., Mol. Cell. Endocrinol. 89:141-151(1992).

