



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

FPRL2 Antibody [#G104]

Catalog Number: G104

Concentration: 1mg/ml

Swiss-Prot No.: P25089

Other Names: FMLP-related receptor II; FMLP-R-II;
Formylpeptide receptor-like 2; FPRL2; FPRH1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human FPRL2.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FPRL2 antibody detects endogenous levels of total FPRL2 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:40000

References:

Bao L., Genomics 13:437-440(1992).

Durstin M., Biochem. Biophys. Res. Commun. 201:174-179(1994).

Kopatz S.A., Submitted (MAR-2003) to the EMBL/GenBank/DDBJ databases.

