



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** +86 519 8805 0026

## FZD1 Antibody [#G105]

**Catalog Number:** G105

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9UP38

**Other Names:** frizzled 1; frizzled-1; Fz-1; FzE1; Wnt receptor

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human FZD1.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** FZD1 antibody detects endogenous levels of total FZD1 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:5000

**References:**

Gazit A., *Oncogene* 18:5959-5966(1999).

Sagara N., *Biochem. Biophys. Res. Commun.* 252:117-122(1998).

Hillier L.W., *Nature* 424:157-164(2003).

