



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

FZD3 Antibody [#G108]

Catalog Number: G108

Concentration: 1mg/ml

Swiss-Prot No.: Q9NPG1

Other Names: Frizzled-3 [Precursor]; Fz-3; hFz3; FZD3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human FZD3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Specificity: FZD3 antibody detects endogenous levels of total FZD3 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 IF: 1:100~1:500


ELISA: 1:1000

References:


Kirikoshi H., Biochem. Biophys. Res. Commun. 271:8-14(2000).

Sala C.F., Biochem. Biophys. Res. Commun. 273:27-34(2000).

Hung B.S., J. Invest. Dermatol. 116:940-946(2001).



Western blot analysis of extracts from RAW264.7 cells and A549 cells, using FZD3 antibody.



Immunofluorescence analysis of COS-7 cells, using FZD3 antibody.