



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

FZD4 Antibody [#G109]

Catalog Number: G109

Concentration: 1mg/ml

Swiss-Prot No.: Q9ULV1

Other Names: Frizzled-4; Fz-4; hFz4; FzE4; CD344 antigen; FZD4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human FZD4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Specificity: FZD4 antibody detects endogenous levels of total FZD4 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%);
Rat (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:10000

References:

Kirikoshi H., Biochem. Biophys. Res. Commun. 264:955-961(1999).

Kopatz S.A., Submitted (NOV-2003) to the EMBL/GenBank/DDBJ databases.

Ota T., Nat. Genet. 36:40-45(2004).

