



Web: [www.anbobio.com](http://www.anbobio.com)

[www.dijibio.com](http://www.dijibio.com)

Order: [order@dijibio.com](mailto:order@dijibio.com)

Tel: (86) 519 8805 0026

## ZNF75 Antibody [#C10689]

**Catalog Number:** C10689

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P51815

**Other Names:** Zinc finger protein 75, Zinc finger protein 82, ZNF75, ZNF82

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human ZNF75.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** ZNF75 antibody detects endogenous levels of total ZNF75 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Villa A., Genomics 18:223-229(1993).

Villa A., Genomics 13:1231-1236(1992).

Marino M., Submitted (AUG-1994) to the EMBL/GenBank/DDBJ databases.

