



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AF4 Antibody [#C10690]

Catalog Number: C10690

Concentration: 1mg/ml

Swiss-Prot No.: P51825

Other Names: AF4, AFF1, FEL protein, MLLT2,
Proto-oncogene AF4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human AF4.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AF4 antibody detects endogenous levels of total AF4 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:20000

References:

Nakamura T., Proc. Natl. Acad. Sci. U.S.A. 90:4631-4635(1993).

Morrissey J., Blood 81:1124-1131(1993).

Gu Y., Cell 71:701-708(1992).



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using AF4 antibody.