



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MAT1 Antibody [#C10691]

Catalog Number: C10691

Concentration: 1mg/ml

Swiss-Prot No.: P51948

Other Names: AMS1, AdoMet synthetase, EC 2.5.1.6, MAT-I/III, MATA1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MAT1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MAT1 antibody detects endogenous levels of total MAT1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Tassan J.-P., EMBO J. 14:5608-5617(1995).

Yee A., Cancer Res. 55:6058-6062(1995).

Tirode F., Mol. Cell 3:87-95(1999).

