



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

NCBP2 Antibody [#C10694]

Catalog Number: C10694

Concentration: 1mg/ml

Swiss-Prot No.: P52298

Other Names: 20 kDa nuclear cap binding protein, CBP20, NCBP 20 kDa subunit, NCBP-interacting protein 1, NCBP2, NIP1, Nuclear cap binding protein subunit 2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

NCBP2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: NCBP2 antibody detects endogenous levels of total NCBP2 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Izaurralde E., Nature 376:709-712(1995).

Lejeune F., EMBO J. 21:3536-3545(2002).

Matsuoka S., Science 316:1160-1166(2007).

