



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

DGKE Antibody [#C10697]

Catalog Number: C10697

Concentration: 1mg/ml

Swiss-Prot No.: P52429

Other Names: Diacylglycerol kinase epsilon, EC 2.7.1.107, Diglyceride kinase epsilon, DGK-epsilon, DAG kinase epsilon, DGKE, DAGK5

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

DGKE.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: DGKE antibody detects endogenous levels of total DGKE protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Tang W., J. Biol. Chem. 271:10237-10241(1996).

Tang W., Gene 239:185-192(1999).

Sjoebloom T., Science 314:268-274(2006).

