



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

GCVK_HHV6Z Antibody [#C10698]

Catalog Number: C10698

Concentration: 1mg/ml

Swiss-Prot No.: P52446

Other Names: Probable ganciclovir kinase, EC=2.7.1.-

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GCVK_HHV6Z.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GCVK_HHV6Z antibody detects endogenous levels of total GCVK_HHV6Z protein.


Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Dominguez G., J. Virol. 73:8040-8052(1999) [PubMed: 10482553].

Lindquister G.J., Arch. Virol. 141:367-379(1996) [PubMed: 8634027].



For Research Use Only