



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MAP2K6 Antibody [#C10699]

Catalog Number: C10699

Concentration: 1mg/ml

Swiss-Prot No.: P52564

Other Names: Dual specificity mitogen-activated protein kinase kinase 6, EC 2.7.12.2, MAP kinase kinase 6, MAP2K6, MAPK/ERK kinase 6, MAPKK 6, MEK6, MP2K6, Mitogen-activated protein kinase kinase 6, PRKMK6, SAPKK3, kinase MKK6

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from N-terminal of human MAP2K6.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MAP2K6 antibody detects endogenous levels of total MAP2K6 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Raingaud J., Mol. Cell. Biol. 16:1247-1255(1996).
Stein B., J. Biol. Chem. 271:11427-11433(1996).
Han J., J. Biol. Chem. 271:2886-2891(1996).

