



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

hnRNP F Antibody [#C10700]

Catalog Number: C10700

Concentration: 1mg/ml

Swiss-Prot No.: P52597

Other Names: HNRPF, ROF, heterogeneous nuclear ribonucleoprotein F

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human hnRNP F.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: hnRNP F antibody detects endogenous levels of total hnRNP F protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Matunis M.J., Nucleic Acids Res. 22:1059-1067(1994).

Honore B., J. Biol. Chem. 270:28780-28789(1995).

McDonald H., Genomics 13:344-348(1992).

