



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ZNF134 Antibody [#C10705]

Catalog Number: C10705

Concentration: 1mg/ml

Swiss-Prot No.: P52741

Other Names: Zinc finger protein 134, ZNF134

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

ZNF134.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ZNF134 antibody detects endogenous levels of total ZNF134 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Tommerup N., Genomics 27:259-264(1995).

