



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## Hexokinase-3 Antibody [#C10708]

**Catalog Number:** C10708

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P52790

**Other Names:** Hexokinase-3, EC 2.7.1.1, Hexokinase type III, HK III, HK3

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human Hexokinase-3.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** Hexokinase-3 antibody detects endogenous levels of total Hexokinase-3 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:5000

**References:**

Furuta H., Genomics 36:206-209(1996).

Palma F., Mol. Cell. Biochem. 155:23-29(1996).

Sjoeblom T., Science 314:268-274(2006).

