



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

EFNA3 Antibody [#C10709]

Catalog Number: C10709

Concentration: 1mg/ml

Swiss-Prot No.: P52797

Other Names: Ephrin-A3 [Precursor], EPH-related receptor tyrosine kinase ligand 3, LERK-3, EHK1 ligand, EHK1-L, EFL-2, EFNA3, EFL2, EPLG3, LERK3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human

EFNA3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: EFNA3 antibody detects endogenous levels of total EFNA3 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Kozlosky C.J., *Oncogene* 10:299-306(1995).

Davis S., *Science* 266:816-819(1994).

Gregory S.G., *Nature* 441:315-321(2006).

