



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## LIMK1 Antibody [#C10719]

**Catalog Number:** C10719

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P53667

**Other Names:** EC 2.7.11.1, KIZ-1, LIM domain kinase 1, LIMK, LIMK-1, kinase LIMK1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human LIMK1.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** LIMK1 antibody detects endogenous levels of total LIMK1 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Mizuno K., *Oncogene* 9:1605-1612(1994).

Frangiskakis J.M., *Cell* 86:59-69(1996).

Osborne L.R., *Genomics* 36:328-336(1996).

