



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## MAPK10 Antibody [#C10722]

**Catalog Number:** C10722

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P53779

**Other Names:** EC 2.7.11.24, JNK3-alpha-2, JNK3A, MAP kinase p49 3F12, MAPK10, MK10, Mitogen-activated protein kinase 10, PRKM10, SAPK-beta, SERK2, Stress-activated protein kinase JNK3, c-Jun N-terminal kinase 3, kinase JNK3, p54- beta

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against

synthesized peptide derived from internal of human MAPK10.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** MAPK10 antibody detects endogenous levels of total MAPK10 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Mohit A.A., Neuron 14:67-78(1995).

Lisnock J., Biochemistry 39:3141-3148(2000).

Jagadish N., Biochem. J. 389:73-82(2005).

