



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

PTTG Antibody [#C10723]

Catalog Number: C10723

Concentration: 1mg/ml

Swiss-Prot No.: P53801

Other Names: C21orf1, PBF, PTTG, PTTG-binding factor, pituitary tumor-transforming gene 1 protein-interacting protein, pituitary tumor-transforming gene protein binding factor

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human

PTTG.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: PTTG antibody detects endogenous levels of total PTTG protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 IHC: 1:50~1:100


ELISA: 1:10000

References:

Yaspo M.-L., Genomics 49:133-136(1998).


Chien W., J. Biol. Chem. 275:19422-19427(2000).

Hattori M., Nature 405:311-319(2000).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human brain tissue using PTTG antibody.



Western blot analysis of extracts from RAW264.7 cells, using PTTG antibody.