



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## EPHB4 Antibody [#C10732]

**Catalog Number:** C10732

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P54760

**Other Names:** EC 2.7.1.112, EC 2.7.10.1, EPB4, Ephrin type-B receptor 4 precursor, HTK, Tyrosine-protein kinase receptor HTK

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

EPHB4.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** EPHB4 antibody detects endogenous levels of total EPHB4 protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Bennett B.D., J. Biol. Chem. 269:14211-14218(1994).

Wilson M.D., Nucleic Acids Res. 29:1352-1365(2001).

Zhang Y., Mol. Cell. Proteomics 4:1240-1250(2005).

