



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AF10 Antibody [#C10737]

Catalog Number: C10737

Concentration: 1mg/ml

Swiss-Prot No.: P55197

Other Names: Protein AF-10, MLLT10, AF10

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human AF10.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AF10 antibody detects endogenous levels of total AF10 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Chaplin T., Blood 85:1435-1441(1995).

Deloukas P., Nature 429:375-381(2004).

Matsuoka S., Science 316:1160-1166(2007).

