



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ELL Antibody [#C10738]

Catalog Number: C10738

Concentration: 1mg/ml

Swiss-Prot No.: P55199

Other Names: RNA polymerase II elongation factor ELL,
RNA polymerase II elongation factor ELL, ELL, C19orf17

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ELL.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ELL antibody detects endogenous levels of total ELL protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Thirman M.J., Proc. Natl. Acad. Sci. U.S.A. 91:12110-12114(1994).

Simone F., Blood 98:201-209(2001).

Kong S.E., Proc. Natl. Acad. Sci. U.S.A. 102:10094-10098(2005).

