



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ADK Antibody [#C10739]

Catalog Number: C10739

Concentration: 1mg/ml

Swiss-Prot No.: P55263

Other Names: ADK2, ATP-AMP transphosphorylase, EC 2.7.4.3, KAD2, adenylate kinase isoenzyme 2, mitochondrial

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human ADK.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ADK antibody detects endogenous levels of total ADK protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Spychala J., Proc. Natl. Acad. Sci. U.S.A. 93:1232-1237(1996).

Singh B., Eur. J. Biochem. 241:564-571(1996).

McNally T., Biochem. Biophys. Res. Commun. 231:645-650(1997).

