



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## P19 INK4d Antibody [#C10741]

**Catalog Number:** C10741

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P55273

**Other Names:** Cyclin-dependent kinase 4 inhibitor D, p19-INK4d, CDKN2D

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human P19 INK4d.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** P19 INK4d antibody detects endogenous levels of total P19 INK4d protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%)


**Applications:** IHC: 1:50~1:100 ELISA: 1:40000

**References:**

Okuda T., Genomics 29:623-630(1995).

Chan F.K.M., Mol. Cell. Biol. 15:2682-2688(1995).

Guan K.L., Mol. Biol. Cell 7:57-70(1996).



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using P19 INK4d antibody.