



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

hnRNP H Antibody [#C10747]

Catalog Number: C10747

Concentration: 1mg/ml

Swiss-Prot No.: P55795

Other Names: FTP-3, FTP3, HNRH2, HNRPH2,
Heterogeneous nuclear ribonucleoprotein H', ROH2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human hnRNP H.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: hnRNP H antibody detects endogenous levels of total hnRNP H protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Vorechovsky I., Genomics 21:517-524(1994).

Ross M.T., Nature 434:325-337(2005).

Olsen J.V., Cell 127:635-648(2006).



For Research Use Only