



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

DLX5 Antibody [#C10748]

Catalog Number: C10748

Concentration: 1mg/ml

Swiss-Prot No.: P56178

Other Names: Homeobox protein DLX-5, DLX5

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human DLX5.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: DLX5 antibody detects endogenous levels of total DLX5 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:40000

References:

Ota T., Nat. Genet. 36:40-45(2004).


Hillier L.W., Nature 424:157-164(2003).

Willis D.M., J. Biol. Chem. 277:37280-37291(2002).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using DLX5 antibody.



Western blot analysis of extracts from HeLa cells, 293 cells and COLO205 cells, using DLX5 antibody.