



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CtBP2 Antibody [#C10751]

Catalog Number: C10751

Concentration: 1mg/ml

Swiss-Prot No.: P56545

Other Names: C-terminal-binding protein 2, CtBP2, CTBP2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CtBP2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CtBP2 antibody detects endogenous levels of total CtBP2 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Katsanis N., Genomics 47:294-299(1998).

Schmitz F., Neuron 28:857-872(2000).

Deloukas P., Nature 429:375-381(2004).

