



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

PER3 Antibody [#C10752]

Catalog Number: C10752

Concentration: 1mg/ml

Swiss-Prot No.: P56645

Other Names: Period circadian protein homolog 3, Circadian clock protein PERIOD 3, Hper3, Cell growth-inhibiting gene 13 protein, PER3, GIG13

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human PER3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: PER3 antibody detects endogenous levels of total PER3 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Gregory S.G., Nature 441:315-321(2006).

Ebisawa T., EMBO Rep. 2:342-346(2001).

Johansson C., Neuropsychopharmacology 28:734-739(2003).

