



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

SOX8/9/17/18 Antibody [#C10757]

Catalog Number: C10757

Concentration: 1mg/ml

Swiss-Prot No.: P57073/P48436/Q9H6I2/P35713

Other Names: Transcription factor SOX-8/9/17/18, SOX8/9/17/18

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human SOX8/9/17/18.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: SOX8/9/17/18 antibody detects endogenous levels of total SOX8/9/17/18 protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:10000

References:

Daniels R.J., Hum. Mol. Genet. 10:339-352(2001).

Martin J., Nature 432:988-994(2004).

The MGC Project Team, Genome Res. 14:2121-2127(2004).




Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human brain tissue using SOX8/9/17/18 antibody.



serum

-

+

-

Peptide

-

-

+

Western blot analysis of extracts from COS7 cells, untreated or treated with serum (10%, 30mins), using SOX8/9/17/18 antibody.