



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

SCAND1 Antibody [#C10759]

Catalog Number: C10759

Concentration: 1mg/ml

Swiss-Prot No.: P57086

Other Names: SCAN domain-containing protein 1, SCAND1, SDP1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human SCAND1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: SCAND1 antibody detects endogenous levels of total SCAND1 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:5000

References:


Schumacher C., J. Biol. Chem. 275:17173-17179(2000).

Deloukas P., Nature 414:865-871(2001).

The MGC Project Team, Genome Res. 14:2121-2127(2004).



Western blot analysis of extracts from Jurkat cells, using SCAND1 antibody.



Immunofluorescence analysis of HepG2 cells, using SCAND1 antibody.