



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TMPRSS3 Antibody [#C10761]

Catalog Number: C10761

Concentration: 1mg/ml

Swiss-Prot No.: P57727

Other Names: Transmembrane protease, serine 3, EC 3.4.21.-, Serine protease TADG-12, Tumor-associated differentially-expressed gene 12 protein, TMPRSS3, ECHOS1, TADG12

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human TMPRSS3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity: TMPRSS3 antibody detects endogenous levels of total TMPRSS3 protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:10000

References:

Underwood L.J., Biochim. Biophys. Acta 1502:337-350(2000).
Scott H.S., Nat. Genet. 27:59-63(2001).
Clark H.F., Genome Res. 13:2265-2270(2003).



Western blot analysis of extracts from HUVEC cells, using TMPRSS3 antibody.



Immunofluorescence analysis of HeLa cells, using TMPRSS3 antibody.