



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

LMO4 Antibody [#C10773]

Catalog Number: C10773

Concentration: 1mg/ml

Swiss-Prot No.: P61968

Other Names: LIM domain only 4, LIM domain transcription factor LMO4, LMO-4, breast tumor autoantigen

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human LMO4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: LMO4 antibody detects endogenous levels of total LMO4 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Racevskis J., Biochim. Biophys. Acta 1445:148-153(1999).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

