



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

BLCAP Antibody [#C10782]

Catalog Number: C10782

Concentration: 1mg/ml

Swiss-Prot No.: P62952

Other Names: Bladder cancer-associated protein,
Bladder cancer 10 kDa protein, Bc10, BLCAP, BC10

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human BLCAP.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Specificity: BLCAP antibody detects endogenous levels of total BLCAP protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%);
Rat (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:40000

References:

Olsen S.R., Mol. Endocrinol. 5:1246-1256(1991).

Knoell R., Submitted (FEB-2000) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team, Genome Res. 14:2121-2127(2004).

