



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## IRX1 Antibody [#C10795]

**Catalog Number:** C10795

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P78414

**Other Names:** ICSBP, ICSBP1, IRF-8, interferon consensus sequence-binding, interferon regulatory factor 8

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human IRX1.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** IRX1 antibody detects endogenous levels of total IRX1 protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:40000

**References:**

Hansen L., Submitted (JUL-2003) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Lewis M.T., Submitted (MAR-1997) to the EMBL/GenBank/DDBJ databases.

