



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

IRX3 Antibody [#C10796]

Catalog Number: C10796

Concentration: 1mg/ml

Swiss-Prot No.: P78415

Other Names: iroquois homeobox 3; IRX-1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human IRX3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: IRX3 antibody detects endogenous levels of total IRX3 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:20000

References:

Hansen L., Submitted (JUL-2003) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Lewis M.T., Submitted (MAR-1997) to the EMBL/GenBank/DDBJ databases.



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using IRX3 antibody.