



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ADARB1 Antibody [#C10800]

Catalog Number: C10800

Concentration: 1mg/ml

Swiss-Prot No.: P78563

Other Names: EC 3.5.-.-, dsRNA adenosine deaminase, RNA-editing deaminase 1, RNA-editing enzyme 1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ADARB1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ADARB1 antibody detects endogenous levels of total ADARB1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Mittaz L., Genomics 41:210-217(1997).

Lai F., Mol. Cell. Biol. 17:2413-2424(1997).

Villard L., Somat. Cell Mol. Genet. 23:135-145(1997).

