



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

FOXK1 Antibody [#C10811]

Catalog Number: C10811

Concentration: 1mg/ml

Swiss-Prot No.: P85037

Other Names: FHX, FOXJ2 forkhead factor, fork head homologous X

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human FOXK1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FOXK1 antibody detects endogenous levels of total FOXK1 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Ota T., Nat. Genet. 36:40-45(2004).

Hillier L.W., Nature 424:157-164(2003).

Katoh M., Int. J. Mol. Med. 14:127-132(2004).

