



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MEF2B Antibody [#C10843]

Catalog Number: C10843

Concentration: 1mg/ml

Swiss-Prot No.: Q02080

Other Names: MDS1-EVI1, PRDM3, myelodysplasia syndrome 1, myelodysplasia syndrome-associated 1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MEF2B.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MEF2B antibody detects endogenous levels of total MEF2B protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Yu Y.-T., Genes Dev. 6:1783-1798(1992).

Pollock R., Genes Dev. 5:2327-2341(1991).

Grimwood J., Nature 428:529-535(2004).

