



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

HEN1/2 Antibody [#C10854]

Catalog Number: C10854

Concentration: 1mg/ml

Swiss-Prot No.: Q02575/Q02577

Other Names: Helix-loop-helix protein 1, HEN1, Nescient helix loop helix 1, NSCL-1, NHLH1, HEN1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human HEN1/2.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: HEN1/2 antibody detects endogenous levels of total HEN1/2 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Brown L., Proc. Natl. Acad. Sci. U.S.A. 89:8492-8496(1992).

Lipkowitz S., J. Biol. Chem. 267:21065-21071(1992).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

