



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

OCT6 Antibody [#C10863]

Catalog Number: C10863

Concentration: 1mg/ml

Swiss-Prot No.: Q03052

Other Names: POU domain, class 3, transcription factor 1, Octamer-binding transcription factor 6, Oct-6, POU domain transcription factor SCIP, POU3F1, OCT6, OTF6

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

OCT6.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OCT6 antibody detects endogenous levels of total OCT6 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Faus I., Mol. Cell. Biol. 14:3263-3275(1994).

Tobler A., Nucleic Acids Res. 21:1043-1043(1993).

