



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CREM Antibody [#C10864]

Catalog Number: C10864

Concentration: 1mg/ml

Swiss-Prot No.: Q03060

Other Names: CREM-alpha, cAMP-responsive element modulator, cAMP-responsive element modulator, alpha isoform CREM, CREM-beta, cAMP-responsive element modulator, cAMP-responsive element modulator, beta isoform CREM-BCEFGIalpha, Tau

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from internal of human CREM.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CREM antibody detects endogenous levels of total CREM protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Masquillier D., Cell Growth Differ. 4:931-937(1993).

Meyer T.E., Nucleic Acids Res. 20:6106-6106(1992).

Fujimoto T., J. Biochem. 115:298-303(1994).

