



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CEBPZ Antibody [#C10870]

Catalog Number: C10870

Concentration: 1mg/ml

Swiss-Prot No.: Q03701

Other Names: CBF, CBF2, CCAAT-binding factor,
CCAAT-box-binding transcription factor

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CEBPZ.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CEBPZ antibody detects endogenous levels of total CEBPZ protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Lum L., Mol. Cell. Biol. 10:6709-6717(1990).

Hillier L.W., Nature 434:724-731(2005).

Olsen J.V., Cell 127:635-648(2006).

