



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

HSF2 Antibody [#C10871]

Catalog Number: C10871

Concentration: 1mg/ml

Swiss-Prot No.: Q03933

Other Names: Heat shock factor protein 2, HSF 2, Heat shock transcription factor 2, HSTF 2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human HSF2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: HSF2 antibody detects endogenous levels of total HSF2 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Schuetz T.J., Proc. Natl. Acad. Sci. U.S.A. 88:6911-6915(1991).

Mungall A.J., Nature 425:805-811(2003).

Sheldon L., Genes Dev. 7:1549-1558(1993).

