



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TLE2 Antibody [#C10874]

Catalog Number: C10874

Concentration: 1mg/ml

Swiss-Prot No.: Q04725

Other Names: Transducin-like enhancer protein 2, ESG2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TLE2.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TLE2 antibody detects endogenous levels of total TLE2 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Stifani S., Nat. Genet. 2:119-127(1992).

Grimwood J., Nature 428:529-535(2004).

Grbavec D., Eur. J. Biochem. 258:339-349(1998).

