



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CSTF1 Antibody [#C10882]

Catalog Number: C10882

Concentration: 1mg/ml

Swiss-Prot No.: Q05048

Other Names: CF-1 50 kDa subunit, CSTF 50 kDa subunit, CSTF1, Cleavage stimulation factor, 50 kDa subunit

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human CSTF1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CSTF1 antibody detects endogenous levels of total CSTF1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Takagaki Y., Proc. Natl. Acad. Sci. U.S.A. 89:1403-1407(1992).

Deloukas P., Nature 414:865-871(2001).

Kleiman F.E., Science 285:1576-1579(1999).

