



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

SRY Antibody [#C10883]

Catalog Number: C10883

Concentration: 1mg/ml

Swiss-Prot No.: Q05066

Other Names: Delayed-early protein HRS, HRS, Pre-mRNA splicing factor SRP40, SFR5, SRP40, Splicing factor, arginine/serine-rich 5

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human SRY.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: SRY antibody detects endogenous levels of total SRY protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Sinclair A.H., Nature 346:240-244(1990).

Su H., Am. J. Hum. Genet. 52:24-38(1993).

Behlke M.A., Genomics 17:736-739(1993).

