



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MEF2C Antibody [#C10891]

Catalog Number: C10891

Concentration: 1mg/ml

Swiss-Prot No.: Q06413

Other Names: Myocyte-specific enhancer factor 2C, MEF2C

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MEF2C.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MEF2C antibody detects endogenous levels of total MEF2C protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:40000

References:

Leifer D., Proc. Natl. Acad. Sci. U.S.A. 90:1546-1550(1993).

McDermott J.C., Mol. Cell. Biol. 13:2564-2577(1993).

Wang A.H., Mol. Cell. Biol. 19:7816-7827(1999).



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using MEF2C antibody.