



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TNFC Antibody [#C10895]

Catalog Number: C10895

Concentration: 1mg/ml

Swiss-Prot No.: Q06643

Other Names: Lymphotoxin-beta, LT-beta, Tumor necrosis factor C, TNF-C, Tumor necrosis factor ligand superfamily member 3, LTB, TNFC, TNFSF3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human TNFC.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TNFC antibody detects endogenous levels of total TNFC protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:10000

References:


Browning J.L., Cell 72:847-856(1993).

Warzocha K., Biochem. Biophys. Res. Commun. 238:273-276(1997).

Neville M.J., J. Immunol. 162:4745-4754(1999).



Western blot analysis of extracts from HepG2 cells, using TNFC antibody.



Immunofluorescence analysis of HeLa cells, using TNFC antibody.