



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

NCBP1 Antibody [#C10911]

Catalog Number: C10911

Concentration: 1mg/ml

Swiss-Prot No.: Q09161

Other Names: CBP80, MGC2087, NCBP, nuclear cap binding protein subunit 1, 80kDa

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human NCBP1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: NCBP1 antibody detects endogenous levels of total NCBP1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Izaurralde E., Cell 78:657-668(1994).

Kataoka N., Nucleic Acids Res. 22:3861-3865(1994).

Humphray S.J., Nature 429:369-374(2004).

